


QUICKIE®

Xtender™

Allows you to
go the distance!



Live Without Limits™

The Quickie® Xtender™ provides power assist to a manual wheelchair, allowing users to maintain the look, control, and functionality of a manual wheelchair. The Xtender increases the distance traveled and the efficiency of propulsion without additional physical exertion.



Powered by Yamaha

- Motorized, power assist wheels increase the force applied to the handrims by 1.5x or 3x
- Electronically synchronized units compensate for arm strength imbalance for optimal response and straight propulsion
- Advanced lightweight battery technology includes a microprocessor for increased performance
- Features rollback assistance for improved performance during hill climbs

Increased Access and Independence

- Reduces energy expenditure by as much as 45%, reduces user fatigue, and increases range of mobility
- Provides excellent handrim grip for users with limited hand function
- May decrease risk of long term shoulder and wrist injuries

Ease of Use

- Quick-release wheels allow for easy transitions between manual and assisted propulsion
- No programming required
- Easy to transport and store
- Bright LED light alerts user when unit is on
- Battery level lights indicate remaining battery power



da Vinci
AWARDS

2008 WINNER

SPECIFICATIONS	Quickie Xtender
HCPCS CODE	E0986
Available On	Quickie 2, GT, GTI, and GP Series
Torque Multiplier	1:1.5 and 1:3
Overall Weight	38 lbs.
Available Handrims	Plastic Coated (1:3) Stainless Steel (1:1.5)
Battery Capacity	7 hours
Time to Fully Charge	2.5 to 3 hours
Maximum Range	9.5 miles
Available Tire	Pneumatic
Weight Capacity	265 lbs.