




# Avid Rehab

## VELOCITY

Avid Rehab's next generation wheelchair provides best in class performance and superior clinical functionality. With enhanced lifestyle elements and state-of-the-art technology the Velocity has easily become Avid Rehab's flagship product.



# Features

- Center of gravity tilt 0° - 50° posterior
- Shear reduction recline 90° - 150°
- Tension-adjustable front and rear suspension
- Now available with  **PRECISION COMFORT** |  **Gel Ovations™**
- Optional power seat lift- 12" telescoping tower
- Adjustable seat width and depth (16" - 22"); adjustments can be made in 1" increments
- Seat to floor height 16.5" - 19.5" (includes power seat lift)
- LED lights and unoccupied transit tie downs included
- Powered by  **MK POWERED**



# Specifications

|   |          |
|---|----------|
| Maximum Speed   | 6 mph    |
| Range   | 18mi     |
| Turning Radius  | 26"      |
| Overall Length  | 36"      |
| Overall Width   | 24"      |
| Ground Clearance                                      | 2.5"     |
| Tilt  | 0°-50°   |
| Recline   | 90°-150° |
| Articulating Foot Platform (seat to footplate length) | 13"-16"  |
| Power Seat Lift                                       | 12"      |
| Seat-to-Footplate Height Captain/Pan Seat             | 10"-13"  |

|                            |   |
|----------------------------|---|
| Drive Wheels               | 14" x 3"                                    |
| Casters Front/Rear         | 6" x 2"                                     |
| Weight Capacity            | 300 lbs                                     |
| Weight Including Batteries | 320 - 404 lbs. (depending on configuration) |
| Electronics                | PG VR2 or R-NET                             |
| Seat Width                 | 16"-22" (by 1" incr.)                       |
| Seat Depth                 | 16"-22" (by 1" incr.)                       |
| Back Height (PC Proform)   | 23"   |
| Seat to Floor Height       | 16.5"-19.5"                                 |
| Batteries                  | 22NF  |
| HCPCS                      | K0848/K0849<br>K0856<br>K0861               |

## Colors



**Avid Rehab**

 Avid\_Rehab

 MeirtsUSA

 Merits Health Products, Inc

